

PERQUE pH Hydrion Papers

Information and Instructions

Directions

Please use only as directed. Withdraw and tear off approximately two inches of paper from the roll and immerse in first morning urine collected after six (6) hours of rest. Remove and promptly compare with specimen colors on dispenser to determine corresponding pH. The target pH range is 6.5 to 7.5.

Size

Measure range of 5.5 – 8. One dispenser equals approximately 90 tests.

14 Pidgeon Hill Drive, Suite 180, Sterling, VA 20165
phone: 800-525-7372 • fax: 703-450-2995 • e-mail: perque@perque.com

PERQUE® is researched, uniquely formulated, and exclusively distributed by PERQUE LLC.